



BD AFFIRM (CANDIDA, GARDNERELLA, TRICHOMONAS)

INDICATIONS:

The BD Affirm Microbial Identification Test is a DNA probe test intended for use in the detection and identification of *Candida* species, *Gardnerella vaginalis* and *Trichomonas vaginalis* nucleic acid in vaginal fluid specimens from patients with symptoms of vaginitis/vaginosis.

SPECIMEN:

Vaginal Sample (BD Affirm Microbial Identification Test)

SUPPLIES:

1. Specimen Requisition
2. Specimen Bag with a biohazard label
3. Vaginal sample (BD Affirm Microbial Transport Swab)

COLLECTION:

1. Open the seal and remove the components. Break the transport reagent ampule and dispense into sample collection tube.
2. Insert an un-lubricated speculum into the vaginal to permit visualization of the posterior vaginal fornix.
3. Using the sterile swab, obtain a sample from the posterior vaginal fornix. Twist or roll the swab against the vaginal wall two or three times ensuring the entire circumference of the swab has touched the vaginal wall. Swab the lateral vaginal wall while removing the swab.
4. Immediately place the swab in the sample collection tube.
5. With the swab touching the bottom of the collection tube, grasp the pre-scored handle of the swab just above the top of the tube and bend until the swab breaks. When the swab is fully inserted into the collection tube, the score mark on the swab is approximately 1 cm. above the top of the collection tube. Discard the broken handle into an infectious waste container.
6. Place the cap over the exposed end of the swab and firmly press the cap onto the tube. The cap will "snap" onto the tube when it is properly seated.
7. Label the sample collection tube with the patient name and date of birth, including the time the sample was collected.
8. Place in a biohazard bag along with the completed requisition in the side pocket of the specimen bag.

STORAGE REQUIREMENTS – Room Temperature

STABILITY REQUIREMENTS – 72 hours



BD AFFIRM (CANDIDA, GARDNERELLA, TRICHOMONAS)

Submit all patient information following the procedure for "Completing a Gynecologic Cytology Requisition".

Complete test requisition including last and first name of patient, patient's date of birth and social security number, body site and source of specimen collected. Label specimen container (using the labels provided on the requisition) with patient's first name and last name, and body site/source. The container must have at least two (2) unique identifiers. Examples of unique identifiers: patient name, DOB, unique bar code, etc. Include pertinent clinical information, i.e., previous malignancy, radiation therapy, drugs, etc. Place container in a specimen bag with a biohazard label. Place the requisition in the side pocket of the specimen bag.