



HERPES I & II (CYTOLOGY)

PURPOSE:

To provide the proper collection protocol for Herpes I & II testing.

MATERIALS:

Cervical Sampling:

- CytoBrush & Spatula or Broom in ThinPrep or SurePath Liquid Base Collection Vial

PROCEDURE:

ThinPrep or SurePath Liquid Based Pap Test:

One specimen will be used to provide samples for both the Pap test and/or the HSV test.

QUALITY CONTROL:

- Minimum Volume:
 - 2 mL
- Storage Requirements:
 - Room Temperature
- Stability Requirements:
 - 3 Weeks